

**Goat Anti Human IgG+IgM+IgA-AbFluor 680**

|                            |  |
|----------------------------|--|
| <b>Catalog No :</b>        | RS4006   |
| <b>Reactivity :</b>        | Human  |
| <b>Applications :</b>      | IF;FCM;WB  |
| <b>Target :</b>            | Human IgG+IgM+IgA  |
| <b>Formulation :</b>       | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide  |
| <b>Source :</b>            | Goat   |
| <b>Dilution :</b>          | Suggested starting 1:50-1:1000 dilutions for most fluorescent applications, while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR® Odyssey equipment would be recommended.   |
| <b>Purification :</b>      | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.  |
| <b>Concentration :</b>     | 1mg/mL   |
| <b>Storage Stability :</b> | Stable for one year at -15°C to -25°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezi   |
| <b>Background :</b>        | Immunoway secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also available and are recommended for use with immunoglobulin-rich samples. |
| <b>Sort :</b>              | 6796   |
| <b>No4 :</b>               | 1  |
| <b>Host :</b>              | Goat   |

---

**Conjugate :** AbFluor 680

---

## Products Images

|                 |         |            |
|-----------------|---------|------------|
| Alexa Fluor 350 | 346/442 | Blue       |
| Alexa Fluor 405 | 401/421 | Blue       |
| Alexa Fluor 488 | 496/519 | Green      |
| Alexa Fluor 532 | 532/553 | Yellow     |
| Alexa Fluor 555 | 555/565 | Yellow     |
| Alexa Fluor 568 | 578/603 | Red/Orange |
| Alexa Fluor 594 | 590/617 | Red/Orange |
| Alexa Fluor 633 | 632/647 | Red        |
| Alexa Fluor 647 | 650/665 | Red        |
| Alexa Fluor 660 | 663/690 | Near IR    |
| Alexa Fluor 680 | 679/702 | Near IR    |
| Alexa Fluor 750 | 749/775 | Near IR    |
| Alexa Fluor 790 | 784/814 | Near IR    |

To use the Alexa Fluors with fluorescent imagers, use a spectral line of the blue laser diode for Alexa Fluors 405, a cyan (488 nm) laser for Alexa Fluors 488, a yellow (526 nm) laser for Alexa Fluor 550 or 594, and a red (633 nm) laser for Alexa Fluor 649. The Alexa Fluor 680 and 790 fluors are compatible with laser- and filter-based infrared imaging instruments that emit in the 700 nm, and 800 nm