

Goat Anti Human IgG Fc-AbFluor 660

Catalog No: RS3904

Reactivity: Human

Applications: Elisa;IF;FCM

Target: Human IgG Fc

Formulation: 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1%

BSA, 50% glycerol, 0.02% Sodium Azide

Source: Goat

Dilution: IF (1:200 - 1:1000), FCM (1:100 - 1:1000), Elisa (Use at an assay dependent

concentration)

Purification: The antibody was isolated from antisera by immunoaffinity chromatography

using antigens coupled to agarose beads.

Concentration: 1mg/mL

Storage Stability: Stable for one year at -15°C to -25°C from date of shipment. For maximum

recovery of product, centrifuge the original vial after thawing and prior to removing

the cap. Aliquot to avoid repeated freezi

Background: Immunoway secondary antibodies are available conjugated to enzyme, biotin or

fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also

available and are recommended for use with immunoglobulin-rich samples.

Products Images

1/2



| Alexa Fluor 350 | 346/442 | Blue |
|-----------------|---------|------------|
| Alexa Fluor 405 | 401/421 | Blue |
| Alexa Fluor 488 | 496/519 | Green |
| Alexa Fluor 532 | 532/553 | Yellow |
| Alexa Fluor 555 | 555/565 | Yellow |
| Alexa Fluor 568 | 578/603 | Red/Orange |
| Alexa Fluor 594 | 590/617 | Red/Orange |
| Alexa Fluor 633 | 632/647 | Red |
| Alexa Fluor 647 | 650/665 | Red |
| Alexa Fluor 660 | 663/690 | Near IR |
| Alexa Fluor 680 | 679/702 | Near IR |
| Alexa Fluor 750 | 749/775 | Near IR |
| Alexa Fluor 790 | 784/814 | Near IR |

To use the Alexa Fluors with fluorescent imagers, use a spectral line of the blue laser diode for Alexa Fluors 405, a cyan (488 nm) laser for Alexa Fluors 488, a yellow (526 nm) laser for Alexa Fluor 550 or 594, and a red (633 nm) laser for Alexa Fluor 649. The Alexa Fluor 680 and 790 fluors are compatible with laser- and filter-based infrared imaging instruments that emit in the 700 nm, and 800 nm

2/2